

WEN Tension Release Comb
acc. to Dr. Wen, DBGM



WEN Tension Release Comb

acc. to Dr. Wen, DBGM


Periodontal interventions for recession coverage are made to cover the root surface as complete as possible. If a soft tissue flap shall be translocated coronally, the WEN Tension Release Comb can be used to extend the split flap by combing, to easily achieve a complete coverage. Depending on the condition of the soft tissue, flap extensions of 15 to 20 mm will be possible.



41.865.01

WEN Tension Release Comb
acc. to Dr. Wen, DBGM,
width of working tips 5 and 8 mm

Application



AESTHETIC IS THE RESULT