

ZEPF Bone Scraper II

In oral surgery, bone replacement material and bone regeneration material is often combined with autogenous bone chips in augmentation.

The Bone Scraper is an ideal instrument for easy and fast harvesting of autogenous bone material. It is used to scrape, collect and transplant autogenous bone.

By using the Bone Scraper, instruments like bone filters, trephines, saws and bone mills are no longer necessary. Its blades allow a harvesting of cortical and cancellous bone.

During scraping, the bone is collected in a chamber and at the same time mixed with blood. The bone material can be implanted directly from the previously sterilized Bone Scraper. The cost-efficient disposable blade offers the practitioner an always sharp and economical instrument.

Essential advantages:

- ergonomically perfect
- easy cleaning
- quickly replaceable disposable blade
- simple application
- cost-efficient











| Illustration | Article Description | Order Quantity |
|--------------|------------------------------------------------------------------------------------------|----------------|
| | 47.957.10 | 1 piece |
| | Bone Scraper II straight, Ø 6 mm with sleeve and one blade, ZEPF -Design | |
| | 47.957.15 | 1 piece |
| | Bone-Scraper II angled, Ø 6 mm with sleeve and one blade, ZEPF -Design | |
| | 47.957.20 | 1 set |
| | Bone-Scraper II straight, Ø 6 mm, with sleeve and blade Ø 6 mm, including 3 spare blades | |
| | 47.957.25 | 1 set |
| | Bone-Scraper II angled, Ø 6 mm, with sleeve and blade Ø 6 mm, including 3 spare blades | |
| | Spare parts for Bone Scraper II | |
| | 47.957.30 | 1 set |
| | 3 Spare Blades, round, Ø 6 mm, plus 1 spare sleeve | |
| | 47.957.50 | 1 set |
| | 5 Spare Blades, round, Ø 6 mm, plus 1 spare sleeve | |
| | 47.957.60 | 1 piece |
| | Handle straight, ZEPF -Design | |
| | 47.957.65 | 1 piece |
| | Handle angled, ZEPF -Design | |
| | 47.957.70 | 1 piece |
| | Sleeve for Bone Scraper II | |









_ AESTHETIC IS THE RESULT





......