


Gracey M5 Curettes Set in Z-SHAPE Handle
with ZEPF nanopal coating



Gracey M5 Curettes Set in Z-SHAPE Handle with ZEPF nanopal coating

With the new ZEPF M5 Gracey Curettes we introduce instruments in a new design.

The Z-SHAPE Relax Handle is designed according to the ergonomic needs of a curettage. The new handle as well as the delicate surfaces of the tip are coated with the new ZEPF nanopal coating. The new ZEPF Gracey Curettes have a 0.9 mm thin, flexible and therefore very delicate 1st shaft. The curette offers a perfect adaption to deep roots. With the nanopal coating the surface hardness is increased to 4500 vickers – as never seen before. The blade is highly quenched and does not need to be sharpened

The instrument tips are available in M4 x 0.5 mm, exchangeable.



24.989.550X

Gracey M5 Z-SHAPE Set, nanopal consisting of:



24.221.01GM5OX

GRA 1/2 M5, yellow



24.225.07GM5OX

GRA 7/8 M5, champagne



24.224.11GM5OX

GRA 11/12 M5, signal purple



24.227.13GM5OX

GRA 13/14 M5, cobalt-blue



24.228.04SOX

Sickle Scaler # 204S, yellow-green

Note:

The M5 Curettes with nanopal coating are also available in ZEPF BIONIK Handle.



AESTHETIC IS THE RESULT